

## Quick-Neuron™ Sensory - Maintenance Medium

Catalog Number: SS-MM

### Introduction

Quick-Neuron™ Sensory - Maintenance Medium may be used for the long-term maintenance of human pluripotent stem cell-derived sensory neurons following differentiation as outlined in the Quick-Neuron™ Sensory - mRNA Kit and Human iPSC-derived Neurons user guides. Quick-Neuron™ Sensory differentiated cell cultures display typical neurite outgrowth and express a variety of neuronal markers, such as tubulin beta 3 class III (TUBB3) and variety of sensory neuron markers such as peripherin (PRPH), islet-1 (ISL1), and brain-specific homeobox/POU domain protein 3A (BRN3A/POU4F1). When handled and maintained according to the instructions in this user guide, sensory neurons are viable long-term and are suitable for a variety of characterization and neurotoxicity assays.

**Scale:** The Quick-Neuron™ Sensory - Maintenance Medium provides sufficient medium for 4 wells of a 24-well plate for up to 2 weeks.

**Related Products:** Quick-Neuron™ Sensory - mRNA Kit, Catalog Number: SM-mRNA

### Kit Contents

Upon receipt, store the reagents at the temperatures indicated in the table below. All reagents are shipped on dry ice.

Reagents	Amount	Storage Conditions
Component N	840 µl	-20°C or -80°C
Component S1	20 µl	-20°C or -80°C
Component P	14 µl	-20°C or -80°C

### Required Consumables

Item	Vendor	Catalog Number
DMEM/F12	ThermoFisher	21331020
Neurobasal Medium	ThermoFisher	21103049
Glutamax (100x)	ThermoFisher	35050061
Penicillin-Streptomycin	ThermoFisher	15140122

### Conditions of Use

This product is for research use only. It is not approved for use in humans or for therapeutic or diagnostic use.

### Technical Support

For technical support, please contact us at [cs@elixirgensci.com](mailto:cs@elixirgensci.com) or call +1 (443) 869-5420 (M-F 9am-5pm EST).

## Media Preparation

---

### Medium N(S1)

1. Prepare Medium N(S1) using the reagents listed in the table below.
  - Thaw Component N on ice for 20-30 minutes.
  - Warm all other reagents at room temperature for 20-30 minutes.
2. Store Medium N(S1) for up to 2 weeks at 4°C. The leftover reagents can be discarded or saved for other uses.

Medium N(S1) Reagents	Volume
DMEM/F12	8 ml
Neurobasal Medium	8 ml
200 mM Glutamax (100x)	84 µl
Penicillin-Streptomycin (10000 units/ml; 100x)	168 µl
Component N	520 µl
Component S1	16.8 µl

### Medium N(S1P)

1. Prepare Medium N(S1P) using the reagents listed in the table below.
  - Thaw Component P at room temperature for 20-30 minutes.

Medium N(S1P) Reagents	Volume
Medium N(S1)	7 ml
Component P	3.5 µl

2. Store Medium N(S1P) for up to 2 weeks at 4°C.
  - The leftover Component P can be discarded or saved for another use.

## First Week

---

1. Warm Medium N(S1P) at room temperature for 20-30 minutes.
2. Pipet out most of the old medium, but not completely (i.e., just enough to cover the surface of the well), from each well using a P1000 pipettor and add 800 µl Medium N(S1P) along the wall of the well very slowly.
3. Incubate the cultures at 37°C, 5% CO<sub>2</sub> for 2 days.
4. For subsequent medium changes, pipet out half (400 µl) of the old medium from each well using a P1000 pipettor and add 400 µl Medium N(S1P).
5. Repeat Step 4 every 2-3 days such as on Monday, Wednesday, and Friday for 1 week.

## Second Week

---

1. Warm Medium N(S1) at room temperature for 20-30 minutes.
2. Pipet out most of the old medium, but not completely (i.e., just enough to cover the surface of the well), from each well using a P1000 pipettor and add 800 µl Medium N(S1) along the wall of the well very slowly.
3. Incubate the cultures at 37°C, 5% CO<sub>2</sub> for 2 days.
4. For subsequent medium changes, pipet out half (400 µl) of the old medium from each well using a P1000 pipettor and add 400 µl Medium N(S1).
5. Repeat Step 4 every 2-3 days such as on Monday, Wednesday, and Friday for 1 week.